

FITC-labeled hCD4(Ectodomain)



Product Name	FITC-labeled hCD4(Ectodomain)
Catalog No	IT-001-CD4p-FITC
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of FITC-labeled protein (0.5 mg/ml) in PBS+40% glycerol
Protein	6xHis tagged human CD4 ectodomain (a.a.1-398) (GenBank No. P01730).
Applications	FACS, protein labeling, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.