




gp120(pWITO4160)(HIV-1)

Product Name	gp120(pWITO4160)
Catalog No	IT-001-RB6p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged HIV-1 gp120 (pWITO4160)(Clade B) protein (a.a.32-505) (Genebank No. AY835451).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(pWITO4160) protein

Note: