




## CD4(*rhesus macaque*)

---

Product Name	CD4 ( <i>rhesus macaque</i> )(a.a.1-398)
Catalog No	IT-001-rhCD4
Source	Ectodomain of <i>rhesus macaque</i> CD4 protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged <i>rhesus macaque</i> CD4 protein (a.a.1-398) (Genebank No. NP_001036127).
Applications	Viral entry inhibition, WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: *rhesus* CD4 protein

Note: