

# HA( $\Delta$ TM)(H1N1)(A/Fort Monmouth/1/1947)



---

Product Name	HA( $\Delta$ TM)(H1N1)(A/Fort Monmouth/1/1947)
Catalog No	IT-003-00102 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM)(H1N1)(A/Fort Monmouth/1/1947) protein (a.a.18-530) (Genbank #: AAA67338)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)(H1N1)(A/Fort Monmouth/1/1947)