




HA₁ (H3N2)

Product Name	HA ₁ (H3N2)(A/Wisconsin/67/X-161/2005 (H3N2))
Catalog No	IT-003-0041p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA ₁ (H3N2)(A/Wisconsin/67/X-161/2005 (H3N2)) protein (a.a.17-346) (Genebank accession #ABO37609)
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HA₁ (H3N2) protein