


HA1(A/Victoria/361/2011)(H3N2)



Product Name	HA1(A/Victoria/361/2011)(H3N2)
Catalog No	IT-003-00423p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA1(A/Victoria/361/2011)(H3N2) protein (a.a. 17-346) (see amino acid sequence below)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ

qklpgndstatlclghavpngtiktndqievtnatelvqnssigeicdsphqildgenctlidallgdpqcdgfnkkw
dlfverskaysncypdyvpyaslrlsvassgtlefnnesfnwtgvtqngtssacirrsnnsffsrlnwlthlnfkypalnvtn
pnneqfdklyiwgvhhpgtdkdqiflyaqssgritvstkrsqqavipnigrprirnipsrisiywtivkpgdillinstgnliap
rgyfkirsgkssimrsdapigkcncsecitpngsipndkpfqnvnritygacpryvkqstlklatgmrnvpekqtrg