


HA1 (A/Indonesia/NIHRD11771/2011)(H5N1)



Product Name	HA1(A/Indonesia/NIHRD11771/2011)(H5N1)
Catalog No	IT-003-00522p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(A/Indonesia/NIHRD11771/2011)(H5N1) protein (a.a. 18-344) (see amino acid sequence below)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ

qicigyhannsteqvdtimeknvtvthaqdilekthngklcdldgvkplilkdcsvagwllgnpmcdefinvpewsyive
kanptndlcypgsfndyeelkhllsrinhfekiqiipksswsdheassgvssacpylgspssffrnvvwlikknstypitkky
nntnqedllilwgihhpnngaeqtmlyqnpttyisigtstlnqrlvpkiatrskvngqsgmeffwtilkpndainfesngnfi
apeyaykivkkgdsaimkseleygncntkcqtpmgainssmpfhnihpitigecpkyvksnrlvlatglrnspqresrrkk