

HA₁(H₈N₄)(Influenza A)



Product Name	HA ₁ (H ₈ N ₄)
Catalog No	IT-003-010p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA ₁ (H ₈ N ₄) (A/pintail duck/Alberta/114/1979(H8N4)) protein (a.a.16-343) (Genebank No. ABB87729)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(H₈N₄) protein

Note: