

# HA<sub>1</sub>(H<sub>10</sub>N<sub>7</sub>)



---

Product Name	HA <sub>1</sub> (A/blue-winged teal/Louisiana/Sg-00073/07(H10N7))
Catalog No	IT-003-011p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA <sub>1</sub> (A/blue-winged teal/Louisiana/Sg-00073/07(H10N7)) protein (a.a. 18-341) (Genebank No. ADH96304)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

