



# HA<sub>1</sub>(H<sub>11</sub>N<sub>2</sub>)(Influenza A)

---

Product Name	HA <sub>1</sub> (H <sub>11</sub> N <sub>2</sub> )
Catalog No	IT-003-012p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA <sub>1</sub> (H <sub>11</sub> N <sub>2</sub> ) (A/duck/Yangzhou/906/2002) protein (a.a.17-342) (Genebank No. AAY85533)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

