



HA₁(H₁₂N₅)(Influenza A)

Product Name	HA ₁ (H ₁₂ N ₅)
Catalog No	IT-003-013p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA ₁ (H ₁₂ N ₅) (A/green-winged teal/ALB/199/1991) protein (a.a.17-342) (Genebank No. ABB88110)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(H₁₂N₅) protein

Note: