



HA(Δ TM)(H₁₃N₈)

Product Name	HA(Δ TM)(H ₁₃ N ₈)
Catalog No	IT-003-014 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(H ₁₃ N ₈)(A/black-headed gull/Netherlands/1/00(H13N8)) protein (a.a.17-530) (Genebank No. AAV91212)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

