



NP (H5N1)

Product Name	NP (H5N1)(A/Indonesia/5/2005(H5N1))	
Catalog No	IT-003-023p	
Source	Viral protein purified from 293 cell expressing NP	
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.5% NP40, 500 mM NaCl, and 30 mM imidazole	
Protein	6xHis tagged nucleocapsid protein(NP) (H5N1) (A/Indonesia/5/2005(H5N1)) protein (Genebank accession #ABI36003)	
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.	
Purity	≥ 95%	
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.	
Use Limitation	For research use only, not for use in diagnostic procedures.	
	<p style="text-align: center;">E</p>	E: purified NP (H5N1) protein