




NP(H3N2)(A/Brisbane/10/2007)

Product Name	NP(H3N2)(A/Brisbane/10/2007)
Catalog No	IT-003-024Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS. Reconstitute the protein with 100µl of Millipore water.
Protein	C-terminal 6xHis tagged NP (H3N2)(A/Brisbane/10/2007) protein (Genebank accession No. EU199248)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at -80°C; <u>AVOID REPEATED FREEZE/THAW</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NP (H3N2)(A/Brisbane/10/2007) protein

Note: