




NS1(A/Brisbane/10/2007 (H3N2))

Product Name	NS1(A/Brisbane/10/2007 (H3N2))
Catalog No	IT-003-031Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS with 8M urea. Reconstitute the protein with 100µl of Millipore water
Protein	C-terminal 6xHis tagged NS1(A/Brisbane/10/2007 (H3N2)) protein (Genebank accession no. ACO95275)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>90%
Storage	Store at -80°C; <u>AVOID FREEZING/THAWING</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NS1(A/Brisbane/10/2007 (H3N2)) protein

Note: