

M1 (A/New Caledonia/20/1999 (H1N1))



| | |
|----------------|--|
| Product Name | M1 (A/New Caledonia/20/1999(H1N1)) |
| Catalog No | IT-003-032Ep |
| Source | Recombinant protein purified from E.coli |
| Formulation | Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water. |
| Protein | C-terminal 6xHis tagged matrix protein 1 (A/New Caledonia /20 /1999(H1N1)) (a.a. 1-252) (Genebank Accession No. ABF21304) |
| Applications | WB standard, antibody ELISA, immunogen, etc. |
| Purity | >95% , based on SDS PAGE |
| Storage | Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

