

# M1(H3N2)(A/Wyoming/3/03)



|                |                                                                                                                                                                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | M1(H3N2)(A/Wyoming/3/03)                                                                                                                                       |
| Catalog No     | IT-003-033Ep                                                                                                                                                   |
| Source         | Recombinant protein purified from E.coli                                                                                                                       |
| Formulation    | Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.                                 |
| Protein        | C-terminal 6xHis tagged M1(H3N2)(A/Wyoming/3/03) matrix protein (a.a. 1-252) (Genebank Accession No. ABW81661).                                                |
| Applications   | WB standard, antibody ELISA, immunogen, etc.                                                                                                                   |
| Purity         | >95% , based on SDS PAGE                                                                                                                                       |
| Storage        | Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.                                                                                                   |

