


HA1(B/Texas/6/2011)



Product Name	HA1(B/Texas/6/2011)
Catalog No	IT-003-B12p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(B/Texas/6/2011) protein (a.a. 16-362) (see amino acid sequence below)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ

driectgitsnsphvvktatqgevntgviplttptksyfanlkgtrtrgkclpcdlnctldldvalgrpmcvgtppsakasilhe
vrpvtsgcfpimhdrtkirqlpnllrgyenirlstqnvidaekapggpyrlgtsgscpnatskigffatmawavpkdnyknt
npltvevpyickeeedqitvwgfhsdnktqmknlygdsnpqkftssangvtthyvsqigdfpdqtedgglpqsgrivvdy
mmqkpgktgtivyqrgvllpqkvwcasgrskvikgslpligeadclhekygglnskpyytgehakaigncpiwvktpl
klangtkyrppakllkerg