



# HA( $\Delta$ TM)(B/Florida/4/2006)

Product Name	HA( $\Delta$ TM)(B/Florida/4/2006)
Catalog No	IT-003-B2 $\Delta$ TMp
Source	Purified protein from transfected 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged ectodomain of human influenza B virus (B/Florida/4/2006) HA protein (a.a.17-546)(Genebank accession no. ACA33493).
Applications	Antibody ELISA, WB standard, etc
Purity	> 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)(B/Florida/4/2006)

Note: