




HA(Δ TM)(B/Florida/4/2006)

Product Name	HA(Δ TM)(B/Florida/4/2006)
Catalog No	IT-003-B2p
Source	Purified protein from transfected 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged ectodomain of human influenza B virus (B/Florida/4/2006) HA protein (a.a.17-546)(Genebank accession no. ACA33493).
Applications	Antibody ELISA, WB standard, etc
Purity	> 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(B/Florida/4/2006)

Note: