

HA₁(B/Wisconsin/01/2010)



Product Name	HA ₁ (B/Wisconsin/01/2010)
Catalog No	IT-003-B7p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA ₁ (B/Wisconsin/01/2010) protein(a.a. 15-361)(Genebank #: EPI271545)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(B/Wisconsin/01/2010)