

HA(Δ TM)(A/New York/6868/2009(H1N1))



Product Name	HA(Δ TM)(A/New York/6868/2009(H1N1))
Catalog No	IT-003-SW15 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/New York/6868/2009(H1N1)) protein (a.a.18-529) (Genebank #: ADD74928)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(A/New York/6868/2009(H1N1))