

# HA(Head Region) (A/Mexico/InDRE4114/2009(H1N1))



---

Product Name	HA(Head Region)(A/Mexico/InDRE4114/2009(H1N1))
Catalog No	IT-003-SW1Hp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA (Head Region) (A/Mexico/InDRE4114/2009(H1N1)) protein (a.a.51-302) (Genebank #: ACQ89903)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

250 kD

150 kD

100 kD

75 kD

50 kD

37 kD

25 kD

20 kD

15 kD

10 kD

E:HA(Head Region)  
(A/Mexico/InDRE4114/2009(H1N1))