



## HA<sub>1</sub>(A/Swine/Iowa/15/1930(H<sub>1</sub>N<sub>1</sub>))

---

Product Name	HA <sub>1</sub> (A/swine/Iowa/15/1930(H <sub>1</sub> N <sub>1</sub> ))
Catalog No	IT-003-SW3p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA <sub>1</sub> (A/swine/Iowa/15/1930(H <sub>1</sub> N <sub>1</sub> )) protein (a.a.18-342) (Genebank No. ABV25634)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA<sub>1</sub>(A/swine/Iowa/15/1930(H<sub>1</sub>N<sub>1</sub>))