


Human Galectin-3/LGALS3 Protein



Product Name	Human Galectin-3/LGALS3 protein
Catalog No	IT-300-004p
Source	Recombinant protein purified from E.coli cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in 50 mM Tris-HCl, 150 mM NaCl, 1 mM EDTA, 0.5 mM TCEP, and 50% Glycerol.
Protein	6xHis tagged Human Galectin-3/LGALS3 protein (a.a.1-250) (GenBank No. P17931).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	> 90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Galectin-3/LGALS3

Note: