

Anti-H1(H1N1)(A/New Caledonia/20/99)



Product Name	Anti- H1(H1N1)(A/New Caledonia/20/99)
Catalog No	IT-003-001M30
Source	Mouse monoclonal antibody, clone 5G1F10
Isotype	IgG1
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing HA(ΔTM)(H1N1)(A/New Caledonia/20/99) protein (IT-003-001ΔTMp).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with HA1(H1N1)(A/New Caledonia/20/99), HA(A/Brisbane/59/2007)(H1N1) and (A/Soloman Island/3/2006)(H1N1), not react with other strains as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

