

Anti-H5(H5N1) (A/Vietnam/1203/2004)



Product Name	Anti- H5(H5N1)(A/Vietnam/1203/2004)
Catalog No	IT-003-0051M1
Source	Mouse monoclonal antibody, clone 3C5
Isotype	IgG2b
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing H5(H5N1)(A/Vietnam/1203/2004) protein (IT-003-0051p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts strongly with H5(A/Vietnam/1203/2004) and (A/Indonesia/5/2005)(H5N1) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

