

Human FGFR3/CD333 Protein (alpha(IIIb))



Product Name	Human FGFR3/CD333 Protein (alpha(IIIb))
Catalog No	IT-300-006p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged Human FGFR3/ACH/CD333/CEK2/JTK4 /HSFGFR3EX protein (a.a.23-377) (GenBank No. NP_001156685).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human FGFR3/CD333 Protein (alpha(IIIb))

Note: