


# Human EpCAM/TROP-1/TACSTD1 Protein



Product Name	Human EpCAM/TROP-1/TACSTD1 Protein
Catalog No	IT-300-007p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged Human EpCAM/TROP-1/TACSTD1 protein (a.a.24-265) (GenBank No. NP_002345).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human EpCAM / TROP-1 / TACSTD1 Protein

Note: