

# Human VGF



Product Name	Human VGF
Catalog No	IT-600-002p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged human VGF (a.a.23-615) (GenBank No. P_003369).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 85%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

