



Mouse CD27 / TNFRSF7

Product Name	Human CD27 / TNFRSF7
Catalog No	IT-300-008mp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged Mouse CD27 / TNFRSF7 protein (a.a.22-182) (GenBank No. NP_001028298.1).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Note:

E: Mouse CD27 / TNFRSF7 Protein