

# HA( $\Delta$ TM)(A/Singapore/INFIMH-16-0019/2016)(H3N2)



Product Name	HA( $\Delta$ TM)(A/Singapore/INFIMH-16-0019/2016)(H3N2)
Catalog No	IT-003-00434 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA( $\Delta$ TM)(A/Singapore/INFIMH-16-0019/2016)(H3N2)(a.a.17-529) (GISAID No. EPI1047604).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq 95\%$
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

