


TpN32 (*T. pallidum*/Nichols strain)



Product Name	TpN32 (<i>T. pallidum</i> /Nichols strain)
Catalog No	IT-501-011Ep
Source	Recombinant protein purified from <i>E. coli</i>
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl with 8M Urea.
Protein	6xHis tagged TpN32 (<i>Treponema pallidum</i> /Nichols strain) (a.a. 1-268) (GenBank Accession No. O07950).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: TpN32 (*T. pallidum*/Nichols strain)