



HA1(A/Hong Kong/125/2017(H7N9))

Product Name	HA1(A/Hong Kong/125/2017(H7N9))
Catalog No	IT-003-0076p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/Hong Kong/125/2017(H7N9)) protein (a.a.19-339) (see below).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Note:

E: HA1 (A/Iowa/32/2016)(H1N2) Protein

Protein SEQ:

dkiclgghavsngrtkvntlttergvevvnatetvertnipricskgkrtvdlgqcglgtitgppqcdqflefsadliierregsdvcypgkfv
neearqilresggidketmgftngirtngvtsackrsgssfyamkwllsntdnaafpqmtksykntrkspaiivwgihhsvstaeqt
klygsgnklvtvgssnyqqsfvpspgarpqvnlgsgridfhwilnnpndtvtfsfngafiapdrasflrgksmgisqvgvqdancegdc
yhsggtiisnlpfqnidravgkcpvrykqrslllatgmknvpeipkgr