



Human OLR1 / LOX1 Protein

Product Name	Human OLR1 / LOX1 Protein
Catalog No	IT-300-055p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human OLR1 / LOX1 Protein (aa 61-273; GenBank# P78380) with 6xHis tagged at C-terminus. The protein forms homodimer at 60-70 kDa.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

