

NA(A/New Caledonia/20/1999(H1N1))



Product Name	NA(A/New Caledonia/20/1999(H1N1))
Catalog No	IT-003-101p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged neuraminidase (A/New Caledonia/20/1999(H1N1)) (aa 36-469; GenBank# ABQ10080).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

