


P30C adhesin (Mycoplasma pneumoniae M129)



Product Name	P30C adhesin (Mycoplasma pneumoniae M129)
Catalog No	IT-503-001Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl, 20% glycerol.
Protein	C-terminal 6xHis tagged P30C adhesin (Mycoplasma pneumoniae M129) (a.a. 98-274) (GenBank Accession No. NP_110141).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: P30C adhesin (Mycoplasma pneumoniae M129)