



## NA(B/Colorado/06/2017)

---

|                |  |
|----------------|--|
| Product Name   | NA(B/Colorado/06/2017)   |
| Catalog No     | IT-003-NB1p  |
| Source         | Recombinant protein purified from 293 cell culture   |
| Formulation    | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.  |
| Protein        | 6xHis tagged neuraminidase (B/Colorado/06/2017)(aa 37-466; GenBank# ASK81308).   |
| Applications   | WB standard, antibody ELISA, Functional analysis, etc  |
| Purity         | ≥95%   |
| Storage        | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |



E: NA(B/Colorado/06/2017)

Note: