



HA (Δ TM)(B/Colorado/06/2017)

Product Name	HA(Δ TM)(B/Colorado/06/2017)
Catalog No	IT-003-B21 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(B/Colorado/06/2017) (aa 18-546; GenBank# ASK81305).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq 90\%$
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

