



Human CXCL12 / SDF-1

Product Name	Human CXCL12 / SDF-1
Catalog No	IT-601-001Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in 50mM Tris-HCl and 150mM NaCl.
Protein	C-terminal 6xHis tagged Human CXCL12 / SDF-1 (GenBank # P48061, aa 22-93)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

