


Human Signal-Regulatory Protein α -1/SIRPA/CD172a



Product Name	Human Signal-Regulatory Protein α -1/SIRPA/CD172a
Catalog No	IT-300-077p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Human Signal-Regulatory Protein α -1/SIRPA/CD172a (aa 31–373; GenBank# P78324).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	\geq 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Signal-Regulatory Protein α -1/SIRPA/CD172a

Note: