


# Human KIR2DL4/CD158d/KIR103



Product Name	Human KIR2DL4/CD158d/KIR103
Catalog No	IT-300-080p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Human KIR2DL4/CD158d/KIR103 (aa 22-242; GenBank# Q99706).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥85%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human KIR2DL4/CD158d/KIR103

Note: