

Human Poliovirus Receptor/PVR/CD155



Product Name	Human Poliovirus Receptor/PVR/CD155
Catalog No	IT-300-082p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Human Poliovirus Receptor/PVR/CD155 (aa 21-343; GenBank# P15151).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Poliovirus Receptor/PVR/CD155

Note: