


Human Vascular Cell Adhesion Protein 1/VCAM-1/CD106/L1CAM



Product Name	Human Vascular Cell Adhesion Protein 1/VCAM-1/CD106/L1CAM
Catalog No	IT-300-091p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Human Vascular Cell Adhesion Protein 1/VCAM-1/CD106/L1CAM (aa 25-698; GenBank# P19320).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Vascular Cell Adhesion Protein 1/VCAM-1/CD106/L1CAM

Note: