


Human Insulin-Like Growth Factor-Binding Protein 4/IGFBP-4



Product Name	Human Insulin-Like Growth Factor-Binding Protein 4/IGFBP-4
Catalog No	IT-401-020p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Insulin-Like Growth Factor-Binding Protein 4/IGFBP-4 (aa 22-258; GenBank# P22692) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Insulin-Like Growth Factor-Binding Protein 4/IGFBP-4

Note: