

Human CD44v10/CD44E/CD44R1



Product Name	Human CD44v10/CD44E/CD44R1
Catalog No	IT-300-102p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human CD44v10/CD44E/CD44R1 (aa 21 - 400; GenBank# P16070-10) with 6xHis tagged at C-terminus. Human CD44E is expressed by cells of epithelium and highly expressed by carcinomas.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human CD44v10/CD44E/CD44R1

Note: