

# Anti-Human CD3 $\epsilon$ /CD3E



Product Name	Anti- Human CD3 $\epsilon$ /CD3E
Catalog No	IT-300-CD3E-m1
Source	Mouse monoclonal antibody, clone 2A9
Isotype	IgG2a
Formulation	Each vial contains 100 $\mu$ g purified IgG in 100 $\mu$ l PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human CD3 $\epsilon$ /CD3E.
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with human CD3 complex as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

