


# Human GDNF Family Receptor $\alpha$ -1/GFRA1



Product Name	Human GDNF Family Receptor $\alpha$ -1/GFRA1
Catalog No	IT-401-023p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	Human GDNF Family Receptor $\alpha$ -1/GFRA1 (aa 25-429; GenBank# P56159) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq$ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human GDNF Family Receptor  $\alpha$ -1/GFRA1

Note: