


Human Transferrin receptor protein 1/TfRC/CD71



Product Name	Human Transferrin receptor protein 1/TfRC/CD71
Catalog No	IT-300-118p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Transferrin receptor protein 1/TfRC/CD71 (aa 89-760; GenBank# P02786) with 6xHis tagged at N-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human Transferrin receptor protein 1/TfRC/CD71

Note: