

Human Advanced glycosylation end product-specific receptor/AGER/RAGE



Product Name	Human Advanced glycosylation end product-specific receptor/AGER/RAGE
Catalog No	IT-600-012p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human RAGE (aa 23-342; GenBank# Q15109) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	

